


## DIMENSIONS



## TECHNICAL DATA

**Total weight of the device:**  
130 Kg (battery included)

**Weight of the towing chassis devoid of the removable control handle:**  
62 Kg

**Weight of the detachable control handle without battery:**  
9,6 Kg

**Weight of the battery box:**  
16,5 Kg

**Maximum capacity:**  
400 Kg

**Minimum viable staircase width:**  
85 cm

**Minimum space for manoeuvre on the landing:**  
Refer to the graphs with the minimum dimensions shown in the catalogue

**Load surface tilting:**  
Automatic and Manual (double operating mode)

**Wheelchair loading platform:**  
Standard dimensions width 70 cm, depth 87.5 cm. Possibility of customisation for different needs.

**Movement on surface:**  
Manual, on sliding wheels/Motorised, by means of independent tracks

**Supply voltage:**  
24 V

**Charging voltage:**  
110-230 V

**Tracks speed:**  
4 metres per minute

**Autonomy with 1 charge:**  
about 300 steps

**Usage temperature:**  
from 5 to 40 °C

**Maximum permissible stairs inclination:**  
35°

**Irreversible gear motors by means of electromagnetic brake**

**Advanced electronics with soft start and proportional joystick**

**Independent track drive:**  
possibility of 360-degree rotation with the tracks

**Increased track section compared to traditional track stairlifts for greater grip and safety.**

\* Zonzini Srl has the right to make changes and / or improvements to the products and / or components at any time, even without prior notice.