

LG 2004 BASIC



New base version

LG 2004 Basic maintains the high safety standards of LG 2004 stairway climbers with tracks, without the aid of 360° rotation. The drawbar does not have folding ends



Stairlift with tracks

The range of stairway climbers with tracks also has a BASIC version that does not offer 360° rotation, it still turns effortlessly on landings thanks to additional small rear wheels. The company has therefore succeeded in achieving the result with its innovative techniques and containing costs.




■ LED for balance point



■ Additional wheels

Positioned at the back, they make manoeuvres on landing easier



■ Safety button protection

Metal protection to avoid accidentally pressing the emergency lock button

Technical Specifications

| | |
|--|---|
| Overall weight: | Kg 44 |
| Base equipment weight | Kg 35 |
| Drawbar weight | Kg 9 |
| Maximum overall dimensions | 150x64x92 Cm |
| Motor power | 300 W |
| Battery charger power supply | 100-240 V |
| Maximum transportable weight | Kg 130 |
| Distance covered in back-up | 540 gradini (con passeggero di 85 Kg) 540 steps (with 85 Kg passenger) |
| Duration with indicator in reserve | 80 gradini / steps |
| Minimum stair width measurement | 77 Cm |
| Minimum manoeuvring space on an L-shaped landing | 98x98 Cm |
| Minimum manoeuvring space on a U-shaped landing | 98x170 Cm |
| Maximum inclination | 70% - 35° |
| Operating temperature | 0°C / 40°C |
| Storage temperature with packaging (dry ambient) | -10° ÷ 50° |
| Maximum speed | 15 Gradini per minuto 15 steps per minute |
| Packaging dimension | 116x70x43 cm |