

LG 2004 RAMP

Stairlift with tracks with functional platforms

The series of stairway climbers with tracks is enriched with its version with ramps for multi-purpose and electric wheelchairs. The maximum capacity is 150 Kg. The folding extensions and the footboard adjustments make wheelchair loading and anchoring easier.










■ Adjustable ramps

Adjustments of 3 cm each



■ Depth-adjustable wheelchair coupling

Adjustment min 16 cm, max 24 cm

LG 2004 Ramp technical specifications

Overall weight:	Kg 84
Base equipment weight	Kg 46
Drawbar weight	Kg 38
Maximum overall dimensions:	144x83x99 cm
Motor power	300 W
Battery charger absorbed power	80 W
Battery charger power supply	100-240 V
Maximum transportable weight	Kg 130 - (optional 150 Kg)
Distance covered in back-up	540 steps (with 85 Kg passenger)
Duration with indicator in reserve	80 gradini / steps
Minimum stair width measurement	93 cm
Minimum manoeuvring space on an L-shaped landing	120x100 cm
Minimum manoeuvring space on a U-shaped landing	120x240 cm
Maximum inclination	70% - 35°
Operating temperature	0°C / 40°C
Storage temperature with packaging (dry ambient)	-10° ÷ 50°
Maximum speed	15 steps per minute
Packaging dimension	116x60x45 cm