



ALUMINIUM WHEELCHAIR

An aluminum wheelchair designed for people who have problems with independent movement, and for people with disabilities and mobility limitations. Powered by the user's hand muscles or guided by an attendant. It can be used both indoors and outdoors. Recommended for amputees. High-quality workmanship, durable construction and the use of innovative solutions guarantee that the trolley will serve for years.

PRODUCT CHARACTERISTICS

- aluminium folding frame
- folding backrest
- attendant brakes

- 5 cm high cushion
- user seat belt
- armrests with a tilting function
- parking brake
- aluminum handrims
- quick release - fast dismantling of wheels
- foldable and removable footrests
- adjustable anti-tip wheels
- reg. center of gravity on the rear and front wheels
- rear pneumatic wheels with a diameter of 24 "
- solid rubber front wheels 8 " in diameter
- frame reinforced with a double cross
- load of 120 kg

TECHNICAL DATA

A Overall width: 62-68 cm
B Inner width: 43, 46, 48 cm
C Overall depth: 111 cm
D Seat height from the ground: 51 (56) cm
E Overall height: 95 (92) cm
FxG Seat dimensions: 36, 39, 41 x 40 cm
Weight: 13.8 kg
Maximum capacity: 120 kg
Folded width: 31 cm
Frame colour: silver
Upholstery: black
 +/- 0.5-1 cm

