



Exclusive


## ALUMINIUM WHEELCHAIR

EXCLUSIVE is an aluminum PREMIUM class wheelchair with a modern design, in a beautiful red colour of the frame, which is intended for people who have problems with independent movement. It works perfectly both indoors as well as outdoors. It is extremely light - it weighs only 13.7 kg. As a standard, the wheelchair is equipped with anti-tip wheels which guarantees exceptional

- pocket on the backrest
- belt for the user
- height adjustable armrests
- adjustment of centre of gravity on front wheels
- quick release
- solid rear wheels 24"
- solid front wheels 8"
- anti-tipping wheels
- reflectors on backrest and wheels

### PRODUCT CHARACTERISTICS

- ergonomic handgrips
- attendant brakes
- folding backrest
- detachable and air permeable upholstery



### TECHNICAL DATA

- A Overall width:** 67 cm
- B Inner width:** 46 cm
- C Overall depth:** 106 cm
- D Seat height from the ground:** 50 cm
- E Overall height:** 91 cm
- FxG Seat dimensions:** 45 x 41 cm
- Weight:** 13.7 kg
- Maximum capacity:** 115 kg
- Folded width:** 34 cm
- Frame colour:** red
- m** +/- 0.5-1 cm