




ALUMINIUM WHEELCHAIR

EXPERT is an aluminum wheelchair for the most demanding users who value extraordinary comfort, practicality and solid workmanship. A wheelchair in a beautiful graphite color frame, which is intended for people who have problems with independent movement. Recommended for amputees due to the adjustable center of gravity on the rear and front wheels. Perfect for home and outdoors. As a standard, the wheelchair is equipped with anti-tip wheels which ensures exceptional comfort and safety of use

PRODUCT CHARACTERISTICS

- removable padded upholstery
- belt for the user

- transverse stabilization
- flip-up and detachable armrests
- backrest angle adjustment
- armrests height adjustment
- seat depth adjustment
- centre of gravity adjustment on the front and rear wheels
- wheel lock
- aluminium handrim
- quick-release
- pneumatic rear wheels 24" in diameter
- solid rubber front wheels 8" in diameter
- swing away footrests
- adjustable anti-tipping wheels
- double crossbar
- reflectors on frame
- pocket on the backrest



TECHNICAL DATA

- A Overall width:** 64-68 cm
B Inner width: 43, 46, 48 cm
C Overall depth: 106 cm
D Seat height from the ground: 50 cm
E Overall height: 92 cm
FxG Seat dimensions: 42, 45, 47 x 42 cm
Weight: 17.5 kg
Maximum capacity: 125 kg
Folded width: 39 cm
Frame colour: graphite
 +/- 0.5-1 cm