



## ALUMINIUM WHEELCHAIR

PREMIUM is an aluminum PREMIUM class wheelchair with a modern design, a detachable upholstery and a folded backrest with a pillow - intended for people who have problems with independent movement. It works perfectly both indoors as well as outdoors. It is extremely light - it weighs up to 15.5 kilograms depending on the size of the seat. As a standard, the wheelchair is equipped with anti-tip wheels which ensures exceptional comfort and safety of use.

**Recommended after amputations.**

### PRODUCT CHARACTERISTICS

- attendant brakes
- (folding) backrest with cushion
- seat belt
- detachable upholstery
- flip-up armrests
- wheel locks
- pneumatic rear wheels 22" in diameter
- front wheels 6" in diameter
- plastic handrims
- quick release - fast dismantling of wheels
- swing away and detachable footrests
- adjustable anti-tip wheels
- frame reinforced with double crossbar

### TECHNICAL DATA

- A Overall width:** 61-66 cm  
**B Inner width:** 43, 46, 48 cm  
**C Overall depth:** 110 cm  
**D Seat height from the ground:** 46 cm  
**E Overall height:** 88 cm  
**F x G Seat dimensions:** 42, 45, 47 cm x 41 cm  
**Weight:** 13.5-15.5 kg  
**Maximum capacity:** 100 kg  
**Folded width:** 36 cm  
 \* depending on seat width

**Frame colour:** silver  
**Upholstery colours:** black and blue, black and claret

+/- 0.5 -1 cm

