

# AMBONA

## 4-wheel walker

### PRODUCT CHARACTERISTICS:

- steel frame
- soft armrest height adjustment
- handle height and angle adjustment
- rear wheels with brakes

### USE:

Walking aid which improves stability and helps to maintain upright ambulation.



Height (cm)

VCBK2A\_01: 97,5 - 137

VCBK2A: 108 - 148



Weight (kg)

10



Maximum capacity (kg)\*

VCBK2A\_01: 125

VCBK2A: 125



Width (cm)

59