



## ALUMINIUM 4-WHEELED ROLLATOR YOLA

The rollator is designed for patients with reduced mobility who need stable support and help while walking. It is also recommended to the elderly as an aid in keeping balance, improving walking comfort, and ensuring safety from falling. The use of the rollator affects positively on the return to activity of the people with locomotion and tilting to erect position problems.

### PRODUCT CHARACTERISTICS

- light aluminium frame
- special crossbar construction
- increased compressive strength
- anatomic handgrips
- height adjustable handles
- hand brakes with locking mechanism
- seat with backrest
- shopping bag
- LED lamp
- crutch holder
- wheels with narrow cross-section and with diameter: 7" for back wheels and 8" for front wheels
- folds into very compact unit for transport

### TECHNICAL DATA

**A Overall width:** 64 cm  
**B Overall depth:** 69 cm  
**CxD Seat dimensions:** 46.5 cm x 26 cm  
**E Seat height from the ground:** 53.5 cm  
**F Handles height adjustment:** 79.5-92.5 cm  
**Weight:** 7.6 kg  
**Maximum capacity:** 130 kg  
**Maximum capacity of the bag:** 5 kg

**Frame colours:** titanium, orange, green  
 +/- 0.5 -1 cm

