



## RECIPROCATING FOLDING WALKER

The walker is designed for patients with reduced mobility who need stable support and help while walking. It is also recommended to elderly as an aid in keeping balance, improving walking comfort, and ensuring safety from falling. The use of the walker affects positively on the return to activity of the people with lower limbs diseases.

### PRODUCT CHARACTERISTICS

- double-reinforced aluminium frame
- wide area of support
- stable and long range of walking
- soft and foamed handgrips
- height adjustable legs
- folds into compact unit for storage and transport
- legs with anti-sliding rubber tips

### TECHNICAL DATA

**A Overall width:** 52 cm

**B Inner width:** 47 cm

**C Overall depth:** 48 cm

**D Height adjustment:** 80-96 cm

**Weight:** 2.8 kg

**Maximum capacity:** 136 kg

**Frame colours:**

JMC-C 3211: brown

JMC-C 3212: silver

+/- 0.5 -1 cm

